

1. Open MP\_DatAn 20.ipf
2. Macros > Compile; then Macros > Initiate the MP DatAn Panel

6 Set the cell # according to experiment notes

4 Enter experiment parameters

7 Set to 0

8

Make sure only Post trace is checked.  
Set the trace number in the "Number" box according to experiment notes.  
Set the "Step" and "# of traces" to 1 to begin with.  
Click "Point #".

Once you have found the traces to use for the IOC panel, change "Step" and "# of traces" to suit your needs.

5 Set experiment date and carry over notes from step 4

3 Set raw ephys file path